




Enzymatic Resistance of Corneas Crosslinked Using Riboflavin in Conjunction With Low Energy, High Energy, and Pulsed UVA Irradiation Modes

Sally Hayes
Structural Biophysics Group
Cardiff University, UK



Schematic of three collagen fibrils showing the likely location of riboflavin/ UVA-induced cross-links



- White square: Within molecules (*intramolecular*)
- Green square: Collagen molecule-collagen molecule at fibril surface (*intermolecular*)
- Blue square: Proteoglycan-collagen molecule (fibril surface)
- Red square: Within proteoglycan core proteins
- Purple square: Proteoglycan core protein-proteoglycan core protein

Method


- Porcine eyes are obtained from the abattoir within ~ 4 hrs of death.
- Randomly divided into cross-linked and non-cross-linked treatment groups.
- After treatment, the cornea is removed and an 8mm disk trephined from the centre.



- Corneal disks are placed in 5ml pepsin digest solution (1g of 500 U/mg pepsin from porcine gastric mucosa in 10ml 0.1 M HCL at pH 1.4) and incubated at 23°C.
- Corneal disk diameter is measured daily until the point of complete digestion.
- Some corneal disks may be removed midway through the digestion process to obtain dry weight measurements.




6 porcine eyes/group	Average time for complete digestion to occur (days)
Untreated	11
Dextran 20%	11
Riboflavin (with 20% dextran)	10
Riboflavin (with 20% dextran) + 3 mW UVA for 30 minutes (SCXL 3 mW)	25



11 porcine eyes/ group	Average time for complete digestion (days)	Average dry weight (undigested tissue mass) at day 12 (g)
Untreated	11	0
Riboflavin+ 3 mW UVA 30 mins (SCXL 3 mW)	25	0.0041
Riboflavin + 9 mW UVA for 10 mins (ACXL 9 mW)	25	0.0020
Riboflavin+ 18 mW UVA for 5 mins (ACXL 18 mW)	25	0.0008

P < 0.0001 (comparing SCXL 3mW to ACXL 9mW and ACXL 18mW)
 P < 0.0001 (comparing ACXL 9mW to ACXL 18mW)



SCXL 3mW > ACXL 9mW > ACXL 18mW

Same diameter but lower undigested tissue mass suggests shallower depth of cross-linking after ACXL.

Also see a shallower demarcation line after ACXL (Kymionis et al. 2014).

11 porcine eyes/group



Untreated

Riboflavin+3mW UVA 30 mins (SCXL 5.4J/cm²)


Riboflavin+30mW UVA 3 mins (HCXL 5.4J/cm²)

Riboflavin+30mW UVA 4 mins (HCXL 7.2J/cm²)

Riboflavin+30mW UVA 8 mins pulsed (p-HCXL 7.2J/cm²)



During high intensity cross-linking the use of pulsed UVA increases the enzymatic resistance of the cornea.



Also see a deeper demarcation line after pulsed high intensity cross-linking (Mazotta et al. 2014).

■ Untreated ■ SCXL (5.4J/cm²) ■ HCXL (5.4J/cm²) ■ HCXL (7.2J/cm²) ■ P-HCXL (7.2J/cm²)

Conclusions

- The intensity and depth of cross-linking varies with different protocols.
- High intensity/same energy protocols result in a shallower depth of cross-linking, possibly due to a more rapid oxygen consumption.
- High intensity/high energy protocols result in more cross-linking in the anterior-most stroma but the depth of cross-linking may be shallower.
- Pulsing UVA during high intensity/high energy procedures can increase the enzymatic resistance of the cornea by increasing oxygen availability.
- The amount of cross-linking needed to stop keratoconus progression is not yet known.

Acknowledgements



Prof Keith Meek
Cardiff University



Nada Aldahlawi
Cardiff University



Prof David O'Brart
St Thomas' Hospital, London

